

bond.lign

Bonding system

Indication

of the visio.lign primers



MKZ-Primer

5 ml

REF MKZ02004



MKZ-Primer

5 ml

REF MKZ02004



**MKZ
EM-Aktivator**

4 ml

REF MKZEM004



K-Primer

4 ml


REF APK25003



visio.link

10 ml

REF VLPMMMA10



visio.lign

The Aesthetic and Functional System

bredent group



Indication of the primers



MKZ-Primer

Preparation of the adhesive bond of composites:

- Titanium
- Non-precious metal alloys
- Zirconium oxide
- Oxide ceramics



MKZ-Primer + MKZ EM-Aktivator (Mix 1:1)

Preparation of the adhesive bond of composites:

- Precious metal alloys (Au, Ag, Pt, Pd)
- eco-alloys (low precious metal alloys)



K-Primer

Preparation of the adhesive bond of composites:

- Lithium (di)silicate
- Veneering and press ceramics

Also suitable for silanization of surfaces.



visio.link

Preparation of the adhesive bond of composites:

- PMMA denture resins
- Composite (veneer composites / composite teeth)
- High-performance polymers BioHPP
- PEEK / PEKK / PAEK

Conditioning of the frameworks

Conditioning of metal and zirconium frameworks (CoCr / NPM / titanium / zirconium)

Sandblast metal frameworks at 3 to 4 bar and zirconium frameworks at max. 2 bar with 110 µm grain aluminium oxide.

The framework must not be cleaned with a steam jet after sandblasting.

Any impurities should be removed using alcohol and a clean brush.

Then apply the MKZ-Primer and wait until it evaporates.

Conditioning of precious metal frameworks (palladium-based / silver-based alloy)

Sandblast metal frameworks with 110 µm grain aluminium oxide at a pressure of 2 to 3 bar. The framework must not be cleaned with a steam jet after sandblasting.

Any impurities should be removed using alcohol and a clean brush. Then mix the MKZ-Primer and the MKZ EM-Aktivator at a ratio of 1:1, apply and wait until the mixture evaporates.

Conditioning of oxide ceramic frameworks (zirconium oxide / aluminium oxide / spinel ceramic)

Sandblast the ceramic frameworks with 110 µm aluminium oxide at a maximum pressure of 2 bar or roughen with a diamond grinder. The framework must not be cleaned with a steam jet after sandblasting/grinding! Remove any impurities using alcohol and a clean brush. Then apply the appropriate primer and wait until it evaporates.

Conditioning of plastics (composite / PMMA materials / high-performance polymers such as BioHPP)

Sandblast the plastic / plastic framework with 110 µm grain aluminium oxide at a pressure of 2 to 3 bar. The framework must not be cleaned with a steam jet after sandblasting. Any impurities should be removed using alcohol and a clean brush. Then apply a thin coating of visio.link and cure for 90 seconds in a light polymerisation device (wavelength range 370 nm - 500 nm). The conditioned area should have a silk-matte finish after light curing. This shows that the layer thickness is perfect.

