

neo.lign

Dientes


Carta de formas





choose natural beauty

powered by
visio.lign

Sup. anteriores


F41	 40,6	 10,2 7,6	
B43	 43,5	 10,4 8,0	
C43	 43,3	 10,7 8,3	
A44	 45,2	 10,2 8,5	
F44	 44,6	 11,0 8,4	
I45	 44,6	 11,3 8,2	
H46	 45,8	 10,2 8,5	

Sup. anteriores

S46		 
I47		 
D48		 
M48		 
D49		 
B51		 

Inf. anteriores


T33




32,8




T35




35,1




D38




38,5




D41







41,9




Sup./Inf. anteriores



Sup./Inf. anteriores





COLORES Y FORMAS "VARIEDAD CON SISTEMA"

neo.lign A

Anterior / dientes anteriores



- 13 formas de maxilar superior
- 4 formas de maxilar inferior

neo.lign P

Posterior / dientes posteriores

Forma G



- Diseño multifuncional de las superficies masticatorias para todos los conceptos de oclusión

Forma L



- Diseño lingualizado de las superficies masticatorias para prótesis totales

Disponible en los

G2 **G3** **G4**

pequeño grande

Disponible en los

L1 **L2** **L3** **L4**

pequeño grande



tablas de combinación neo.lign Dientes

Anteriores		Posteriores Sup. / Inf.	
Sup.	Inf.	Forma G	Forma L
F41	T33	G2	L1 / L2
B43	T33 / T35	G2	L2
C43	T33 / T35	G2	L2
A44	T33 / T35	G2 / G3	L2 / L3
F44	T33 / T35	G2 / G3	L2 / L3
I45	T33 / T35	G2 / G3	L2 / L3
H46	T33 / T35	G2 / G3	L2 / L3
S46	T35 / D38	G3	L3
I47	D38	G3	L3
D48	D38	G3 / G4	L3 / L4
M48	D38	G3 / G4	L3 / L4
D49	D38 / D41	G3 / G4	L3 / L4
B51	D41	G4	L4

neo.lign Dientes a novo.lign Carillas Posteriores Sup. / Inf.

neo.lign Forma L	novo.lign Forma W
L1	-
L2	W3
L3	W4
L4	W5

neo.lign Carta de formas

Dientes



choose natural beauty

bre^{dent}
group

0003290E-20210629
Salvo error y modificaciones

