

			Shades (classic A-D shade system)														
novo.lign A			BL3	A1	A2	A3	A3,5	A4	B2	B3	B3	C2	C3	D2	D3	D4	
Upper anterior veneers			REF	BL3	A10	A20	A30	A35	A40	B20	B30	B40	C20	C30	D20	D30	D40
C43	Set of 6 (13, 12, 11, 21, 22, 23)	VOC43	---	5													
C43	Set of 4 (12, 11, 21, 22)	VOC43	---	4													
C43	Set of 2 (13, 23)	VOC43	---	3													
A44	Set of 6 (13, 12, 11, 21, 22, 23)	VOA44	---	5													
A44	Set of 4 (12, 11, 21, 22)	VOA44	---	4													
A44	Set of 2 (13, 23)	VOA44	---	3													
I45	Set of 6 (13, 12, 11, 21, 22, 23)	VOI45	---	5													
I45	Set of 4 (12, 11, 21, 22)	VOI45	---	4													
I45	Set of 2 (13, 23)	VOI45	---	3													
S46	Set of 6 (13, 12, 11, 21, 22, 23)	VOS46	---	5													
S46	Set of 4 (12, 11, 21, 22)	VOS46	---	4													
S46	Set of 2 (13, 23)	VOS46	---	3													
I47	Set of 6 (13, 12, 11, 21, 22, 23)	VOI47	---	5													
I47	Set of 4 (12, 11, 21, 22)	VOI47	---	4													
I47	Set of 2 (13, 23)	VOI47	---	3													
D48	Set of 6 (13, 12, 11, 21, 22, 23)	VOD48	---	5													
D48	Set of 4 (12, 11, 21, 22)	VOD48	---	4													
D48	Set of 2 (13, 23)	VOD48	---	3													
M48	Set of 6 (13, 12, 11, 21, 22, 23)	VOM48	---	5													
M48	Set of 4 (12, 11, 21, 22)	VOM48	---	4													
M48	Set of 2 (13, 23)	VOM48	---	3													
D49	Set of 6 (13, 12, 11, 21, 22, 23)	VOD49	---	5													
D49	Set of 4 (12, 11, 21, 22)	VOD49	---	4													
D49	Set of 2 (13, 23)	VOD49	---	3													
B51	Set of 6 (13, 12, 11, 21, 22, 23)	VOB51	---	5													
B51	Set of 4 (12, 11, 21, 22)	VOB51	---	4													
B51	Set of 2 (13, 23)	VOB51	---	3													
K53	Set of 6 (13, 12, 11, 21, 22, 23)	VOK53	---	5													
K53	Set of 4 (12, 11, 21, 22)	VOK53	---	4													
K53	Set of 2 (13, 23)	VOK53	---	3													
<b>Lower anterior veneers</b>																	
T35	Set of 6 (43, 42, 41, 31, 32, 33)	VUT35	---	5													
T35	Set of 4 (42, 41, 31, 32)	VUT35	---	4													
T35	Set of 2 (43, 33)	VUT35	---	3													
D38	Set of 6 (43, 42, 41, 31, 32, 33)	VUD38	---	5													
D38	Set of 4 (42, 41, 31, 32)	VUD38	---	4													
D38	Set of 2 (43, 33)	VUD38	---	3													
D41	Set of 6 (43, 42, 41, 31, 32, 33)	VUD41	---	5													
D41	Set of 4 (42, 41, 31, 32)	VUD41	---	4													
D41	Set of 2 (43, 33)	VUD41	---	3													

			Shades (classic A-D shade system)														
novo.lign P, multifunctional			BL3	A1	A2	A3	A3,5	A4	B2	B3	B4	C2	C3	D2	D3	D4	
Veneers, posterior G 3			REF	BL3	A10	A20	A30	A35	A40	B20	B30	B40	C20	C30	D20	D30	D40
1G3	1. quadrant, (14, 15, 16, 17)	VO1G3	---														
2G3	2. quadrant, (24, 25, 26, 27)	VO2G3	---														
3G3	3. quadrant, (34, 35, 36, 37)	VU3G3	---														
4G3	4. quadrant, (44, 45, 46, 47)	VU4G3	---														
<b>Veneers, posterior G 4</b>																	
1G4	1. quadrant, (14, 15, 16, 17)	VO1G4	---														
2G4	2. quadrant, (24, 25, 26, 27)	VO2G4	---														
3G4	3. quadrant, (34, 35, 36, 37)	VU3G4	---														
4G4	4. quadrant, (44, 45, 46, 47)	VU4G4	---														
<b>Veneers, posterior W 3</b>																	
1W3	1. quadrant, (14, 15, 16, 17)	VO1W3	---														
2W3	2. quadrant, (24, 25, 26, 27)	VO2W3	---														
3W3	3. quadrant, (34, 35, 36, 37)	VU3W3	---														
4W3	4. quadrant, (44, 45, 46, 47)	VU4W3	---														

Please enter the order quantity.

combo.lign	REF	BL3	A1	A2	A3	A3,5	A4	B2	B3	B4	C2	C3	D2	D3	D4
		BL3	A10	A20	A30	A35	A40	B20	B30	B40	C20	C30	D20	D30	D40
Fixation composite, 8 g	CO2X4 ---														
Mixing cannulas, 10 pcs.	COM KG21 0														

combo.lign gingiva colored	Cont.	REF	Quantity
Dual-hardening composite	8 g	CO2X4GUM	

visio.link	Cont.	REF	Quantity
PMMA and Composite Primer	10 ml	VLP MMA1 0	

visio.sil	Cont.	REF	Quantity
Transparent A - silicone	50 ml	540 0120 0	
Mixing cannulas, visio.sil	12 Pcs.	320 0045 7	
visio.sil fix	50 ml	540 0130 0	
Mixing cannulas, visio.sil	12 Pcs.	320 0045 1	

combo.lign Opaque	Cont.	REF	Quantity
Opaque light, light-curing	4 g	CO1X4OPL	
Opaque medium, light-curing	4 g	CO1X4OPM	
Opaque intensiv, light-curing	4 g	CO1X4OPI	
Opaque GUM, light-curing	4 g	CO1X4OPG	
Opaque catalyst, self-curing	4 g	CO1X4Kat	
Starterkit, all colors, 2 x Kat	6 x 4 g	OLIGNSET1	

beauty setup wax	Cont.	REF	Quantity
beauty setup bright	25 g	430 0030 0	
beauty setup dark	25 g	430 0031 0	

Dispenser	Cont.	REF	Quantity
Dispenser 5 ml 1:1	1 Pc.	320 0044 1	

visio.lign Toolkit	Cont.	REF	Quantity
Composite processing kit	1 Pc.	VLTOOLKIT	

crea.lign accessories	Cont.	REF	Quantity
Handle for turn-type syringe	1 Pc.	320 0044 2	
Application cannulas	10 Pcs.	320 0094 0	
crea.lign Syringe holder	1 Pc.	320 0044 3	

crea.lign	Cont.	REF	Quantity
crea.lign Dentine BL3	5 g	CLFNDBL3	
crea.lign Dentine A1	5 g	CLFNDA10	
crea.lign Dentine A2	5 g	CLFNDA20	
crea.lign Dentine A3	5 g	CLFNDA30	
crea.lign Dentine A3,5	5 g	CLFNDA35	
crea.lign Dentine A4	5 g	CLFNDA40	
crea.lign Dentine B2	5 g	CLFNDB20	
crea.lign Dentine B3	5 g	CLFNDB30	
crea.lign Dentine B4	5 g	CLFNDB40	
crea.lign Dentine C2	5 g	CLFNDC20	
crea.lign Dentine C3	5 g	CLFNDC30	
crea.lign Dentine D2	5 g	CLFNDD20	
crea.lign Dentine D3	5 g	CLFNDD30	
crea.lign Dentine D4	5 g	CLFNDD40	

crea.lign Incisal E2	5 g	CLFN00E2	
crea.lign Incisal E3	5 g	CLFN00E3	
crea.lign Incisal opal	5 g	CLFN00I1	
crea.lign Incisal blue	5 g	CLFN00I2	
crea.lign Incisal rose	5 g	CLFN00I3	
crea.lign Incisal universal	5 g	CLFN00I4	

crea.lign GUM light	5 g	CLFN00G1	
crea.lign GUM rosa	5 g	CLFN00G2	
crea.lign GUM pink	5 g	CLFN00G3	

crea.lign Modifier beige	5 g	CLFN00M1	
crea.lign Modifier oliv	5 g	CLFN00M2	
crea.lign Modifier caramel	5 g	CLFN00M3	
crea.lign Modifier lila	5 g	CLFN00M4	

crea.lign Stain orange	1,4 g	CLFN00S1	
crea.lign Stain brown	1,4 g	CLFN00S2	
crea.lign Modelling Liquid	10 ml	CLFMOD10	

crea.lign Starter kit													REF		Quantity	
10 syringes, 5 g each (at customer's option)													CLIGNSETN			
DBL3	DA1	DA2	DA3	DA3,5	DA4	DB2	DB3	DC3	DD3	E2	E3	I opal	I blue	G light	G rosa	G pink

Please mark 10 shades with X

visio.lign veneering system Starter kit 2	REF	Shade
10 x novo.lign A, 8 x novo.lign P, 1 x combo.lign, 1 x Mixing cannulas, 1 x visio.link, 21 pieces, shade at customer' s option	VLIGNSET2	

visio.lign veneering system Starter kit 3	REF	Shade
such as VLIGNSET2, 2 x crea.lign, 2x Opaque combo.lign, 12 (instead 10) x novo.lign A, 27 pieces, shade at customer' s option	VLIGNSET3	

Sender (Stamp):

Customer No.

Date, Signature

neo.lign A Anterior teeth, upper jaw		REF	Shades (classic A-D shade system)															
			A1	A2	A3	A3,5	A4	B2	B3	B3	C2	C3	D2	D3	D4			
			A10	A20	A30	A35	A40	B20	B30	B40	C20	C30	D20	D30	D40			
C43	Set of 6 (13, 12, 11, 21, 22, 23)	TOC43	---	S														
A44	Set of 6 (13, 12, 11, 21, 22, 23)	TOA44	---	S														
I45	Set of 6 (13, 12, 11, 21, 22, 23)	TOI45	---	S														
S46	Set of 6 (13, 12, 11, 21, 22, 23)	TOS46	---	S														
I47	Set of 6 (13, 12, 11, 21, 22, 23)	TOI47	---	S														
D48	Set of 6 (13, 12, 11, 21, 22, 23)	TOD48	---	S														
M48	Set of 6 (13, 12, 11, 21, 22, 23)	TOM48	---	S														
D49	Set of 6 (13, 12, 11, 21, 22, 23)	TOD49	---	S														
B51	Set of 6 (13, 12, 11, 21, 22, 23)	TOB51	---	S														
<b>Anterior teeth, lower jaw</b>																		
T35	Set of 6 (43, 42, 41, 31, 32, 33)	TUT35	---	S														
D38	Set of 6 (43, 42, 41, 31, 32, 33)	TUD38	---	S														
D41	Set of 6 (43, 42, 41, 31, 32, 33)	TUD41	---	S														

neo.lign P, multifunctional Posterior teeth, size 2		REF	Shades (classic A-D shade system)															
			A1	A2	A3	A3,5	A4	B2	B3	B3	C2	C3	D2	D3	D4			
			A10	A20	A30	A35	A40	B20	B30	B40	C20	C30	D20	D30	D40			
1G2	1. quadrant, (14, 15, 16, 17)	T01G2	---															
2G2	2. quadrant, (24, 25, 26, 27)	T02G2	---															
3G2	3. quadrant, (34, 35, 36, 37)	TU3G2	---															
4G2	4. quadrant, (44, 45, 46, 47)	TU4G2	---															
<b>Posterior teeth, size 3</b>																		
1G3	1. quadrant, (14, 15, 16, 17)	T01G3	---															
2G3	2. quadrant, (24, 25, 26, 27)	T02G3	---															
3G3	3. quadrant, (34, 35, 36, 37)	TU3G3	---															
4G3	4. quadrant, (44, 45, 46, 47)	TU4G3	---															
<b>Posterior teeth, size 4</b>																		
1G4	1. quadrant, (14, 15, 16, 17)	T01G4	---															
2G4	2. quadrant, (24, 25, 26, 27)	T02G4	---															
3G4	3. quadrant, (34, 35, 36, 37)	TU3G4	---															
4G4	4. quadrant, (44, 45, 46, 47)	TU4G4	---															

Design chart neo.lign (full) teeth	REF	Shade
24 pieces (12 x neo.lign A, 12 x neo.lign P), shade at customer's option	VLIGNPRV	

neo.lign Upgrade Kit													REF	Quantity		
24 pieces (12 x neo.lign A, 12 x neo.lign P)													NLUPSET2			
DBL3	DA1	DA2	DA3	DA3,5	DA4	DB2	DB3	DC3	DD3	I rose	I opal	I blue	I univ.	G light	G rosa	G pink

Please enter the order quantity.

Sender (Stamp):

Customer No.

Date, Signature